**B W / A**  
BE WHO YOU ARE

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**ENG** Slight itching is normal as your skin and hair adjust to being covered. If you are experiencing extreme discomfort, itching or pain, remove trans tape as soon as you can, safely. Overstretching as well as improper removal are the most common causes of skin irritation or blistering.

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**Composition:** TransTape is made from a blend of 95% cotton and 5% spandex with a hypoallergenic medical grade acrylic adhesive that is latex free and non-medicated.

**Find more tips and videos for «taping» on our website in the «About» section.**

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Resistant  
Comfortable  
Waterproof

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
→ Apply the tape on the first time. If you peel it off to stick it back on, it will no longer stick.

→ TRANS TAPE is waterproof and sweatproof so it can be worn in the shower, the pool and even in the ocean!

## Caution

### 1 Before applying the tape

Make sure you do not react to the material. Cut a small end and place it on your forearm. Wait 8 hours and check that there is no reaction.

Make sure your skin is clean and dry. Create 2 small squares of soft paper or cotton the size of your nipples. Place them on your nipples before applying the tape and fix them with a short stripe of tape.

### 2 Apply the tape

The tape will stretch a lot, so start 1/3 shorter than you think.

Cut all strips to size, then round the corners with scissors to avoid curling up on the edges and to create a more contoured look.

You don't want it to stretch onto your back, just under your armpit. You should only need two pieces per side, three maximum. Less is better. Always start on the inside of your chest and stretch the piece to under your armpits.

### 3 To remove the tape

Apply our sweet almond oil solution on and all around the tape. Massage the tape and skin to allow the oil to soak through for approximately 1-3 minutes before next step. Start to remove the tape from the inside of your chest, then gradually remove the rest it.

Once removed, massage the area with the remaining oil and apply a moisturizer.

### 4 Care

If your skin is sensitive you may develop redness. Never re-tape on an irritated area, wait until your skin has healed before taping again.